

LINE CATHETER

Epidural Catheter System



www.zenitht.co.kr



ZENITHT

LINE CATHETER


Epidural Catheter System

● Indication

- Post-laminectomy Syndrome
- Chronic back pain
- Spinal Stenosis
- Epidural Adhensiolysis
- Failed Back Surgery Syndrome
- Herniated Nucleus Pulposus

● Features

- The stylet is making the soft catheter less soft and to cut the dura mater more easily
- No need to withdrawn the stylet to inject the medication
- Surgical grade stainless steel
- Flexible atraumatic tip design
- Kink & collapse resistant
- Radiopaque for distinct images and placement acuracy
- Non-reactive with neruolytics : Phenol, alchol, glycerol, etc.
- Tensile spring guided catheter
- Stimulation



Product	LINE CATHETER	
Model	Line Catheter-L	Line Catheter-C
Length	610 mm	305 mm
Diameter	1.2 mm	0.8 mm