

JENITH- Balloon




JENITH-Balloon

● Indication

- Spinal stenosis
- Herniated Nucleus Pulposus
- Chronic Back Pain
- Adhesion & Inflammation
- Post-Laminectomy Syndrome

● Features

- Perform bidirectionally
- Radiopaque & soft aromatic Shaft
- Steerable atraumatic distal tip
- Easy to handle
- Easy to steer precisely
- No need to withdraw the guidedwire to inject the medication



	Overall Length	Catheter Length	Catheter Diameter	Balloon Diameter (Max)
Catheter	485 mm	300 mm	1.7 mm / 2.0 mm	10 mm
Needle	110 mm	85 mm	2.7 mm / 3.0 mm	-