

JENITH-L



ZENITHT



www.zenitht.co.kr


JENITH-L

● Indication

- Spinal stenosis
- Spondylolisthesis and radiculopathy
- Back pain with or without radiculopathy
- Failed back surgery syndrome
- Disk herniation and radiculopathy
- Radiculopathy due to Epidural Fibrosis
- Discogenic back pain and back spasm

● Features

- Perform bidirectionally
- Radiopaque & soft aromatic Shaft
- Steerable atraumatic distal tip
- Easy to handle
- Easy to steer precisely
- No need to withdrawn the stylet to inject the medication



	Overall Length	Catheter Length	Catheter Diameter
Catheter	485 mm	300 mm	1.7 mm / 2.0 mm
Needle	110 mm	85 mm	2.4 mm / 2.7 mm