

# JENITH-SL



ZENITHT



[www.zenitht.co.kr](http://www.zenitht.co.kr)

# JENITH-SL

## ● Indication

- A. Chronic Back Pain
- B. Post-Laminectomy Syndrome
- C. Spinal stenosis
- D. Herniated Nucleus Pulposus

## ● Features

- Lightweight, ergonomic design
- Catheter outer diameter is 3mm
- Soft, atraumatic and radiopaque tip
- Dual working channels
- Dual infusion ports with rubber buffers
- Direct visualization via endoscope
- 2-way steering
- Locking / Unlocking button



	Overall Length	Catheter Length	Catheter O.D	Catheter I.D
Catheter	485 mm	300 mm	3.0 mm	1.3 mm
Compatible	Endoscope		1.2 mm 이하	
	Laser		1.2 mm 이하	