

ZEUS

KYPHO PLASTY Balloon Catheter System

ZEUS-Balloon Specification

NO.	Part NO.	Initial Length (mm)	PSI Rating	Max. Inflation Volume(cc)	Maximum Diameter(mm)	Maximum Length(mm)
1	ZE-BA10	10	300	3	14	16
2	ZE-BA15	15	300	4	17	22
3	ZE-BA20	20	300	6	19	34

ZEUS-KIT Composition

NO.	Description	1 Set	2 Set
1	Cannula	1ea	2ea
2	Bone Drill	1ea	2ea
3	Bone Cement Filler	3ea	6ea
4	Guide Wire	1ea	2ea

ZEUS-Needle (11G, 12G, 13G)

NO.	Part NO.	Description	Size
1	ZE-VP11	VP NEEDLE	11G
2	ZE-VP12		12G
3	ZE-VP13		13G

ZEUS

KYPHO PLASTY Balloon Catheter System

Indication

- VCF(Vertebral compression fracture)due to Osteoporosis
- Osteolytic fracture
- Metastatic bone fractures



ZEUS

KYPHO PLASTY Balloon Catheter System


• ZEUS-Balloon Balloon Catheter 10 / 15 / 20



Special Feature

- Easy to inflate and Variety size


• Balloon Expander Balloon Inflation Device



Special Feature

- Easy to inflate and Fluorescence

• ZEUS-KIT



• ZEUS-Needle 11G, 12G, 13G

